



BTS Provides the Best tools for Protein Purification!

AGAROSE BEADS ACCESSORIES

TECHNICAL SPECIFICATIONS Empty Columns

COLUMN MATERIAL	Polypropylene	
FRIT MATERIAL	Polyethylene	
FRIT PORE SIZE	20 μ m	
CHEMICAL STABILITY	Stable in all commonly used reagents	
PRODUCT NAME	Plastic Columns 50 columns for single-use applications	Plastic Columns XL 50 columns for single-use applications
CAT. No.	C-50	CXL-50
CAPACITY (ml)	12	35