



CERTIFICATE OF ANALYSIS

Catalog No.: BRPM205-50, 205-250


Product Name: Pretained SDS-PAGE standards (Broad range, 2-200kDa)

Description: 8 protein bands. Prestained SDS-PAGE Standard contains seven protein bands: 200kDa, 105kDa, 71kDa, 50kDa, 35kDa, 16kDa, and 2kDa. These proteins provide a set of unambiguous markers for gel electrophoresis. Furthermore, the protein markers can be used in conjunction with western blotting technique to provide a visual display of marker proteins on the transfer membrane.

Recommended

Take out the standards from refrigerator, and thaw. Do not boil. For western blotting, apply approximately 5 μ l of standards per well in mini-gel, or 10 μ l of standards per well in full length-gel. For visual aid, apply approximately 10-15 μ l of standards per well in mini-gel, or 20-30 μ l of standards per well in full length-gel.

Quantity: Sufficient for 50 & 250 loadings mini-gels
(8x10cm) Storage: -20 °C



R X SDS PAGE

Product Size Catalog No. Protein Markers:

50 Loadings BRPM210 Pre stained, Wide Range

250 Loadings BRPM209

Quality authorized by: Date of Analysis: Lot: